



ISO 9001:2000
ISO 13485:2003

Sure Ortho Light Cure Retainer

TECHNICAL DATA

Adhesive bonding to etched enamel	> 30 MPa
Adhesive bonding to Rexillum	> 40 MPa
Water sorption	30 $\mu\text{g}/\text{mm}^3$
Solubility	7 $\mu\text{g}/\text{mm}^3$
Film thickness	15 μm
Dimensional Change on Polymerization	3 %
Working Time (<i>in ambient light and temperature</i>)	2 - 3 min
Setting Time	10 - 20 sec
Compatible with halogen light	Yes
Compatible with Plasma ark light lamp	Yes
Compatible with LED	Yes
Shelf Life	2 Years

||| SURE
orthodontics

Sure Orthodontics
P.O. Box 3362 - 1211 Geneva 3 - Switzerland